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Zdroje:

1. Gilchrest BA, Krutmann J. Skin Aging. Heidelberg: Springer Science & Business Media; 2006. Zdroj: <https://www.sciencedirect.com/science/article/pii/S0190962289702279>
2. Shu YY, Maibach HI. Estrogen and Skin. Am J Clin Dermatol. 2011;12(5):297–311. doi: 10.2165/11589180-000000000-00000 Zdroj: <https://pubmed.ncbi.nlm.nih.gov/21714580/>
3. Farage MA, Miller KW, Maibach HI. Textbook of Aging Skin. Heidelberg: Springer Science & Business Media; 2009. Zdroj: https://www.researchgate.net/publication/239619909_Textbook_of_Aging_Skin
4. Baek, H. J., & Lee, J.-N. (2020). Influence of Value Consciousness on Purchase Behavior and Repurchase Intention for Personal Beauty Devices among Women in Self-Beauty Care. Journal of the Korean Applied Science and Technology, 37(4), 958-968. Zdroj: https://www.researchgate.net/publication/373547984_The_Influence_of_Consumption_Value_Consciousness_on_Beauty_Equipment_Selection_Attributes_and_Purchase_Intention_of_Adults_in_Fujian_Province_of_Southern_China
5. Lim, K. H., Ryu, J. Y., An, I. S., Kim, G. Y., & Lee, J. M. (2009). Anti-wrinkle effects of irregularity low-frequency stimulation. Kor J Aesthet Cosmetol, 7(2), 137-147. Zdroj: [불규칙저주파를 이용한 주름개선효과 임계화.pdf](#)
6. Khakhar, A. Facial EMS Is No Longer a Fringe Therapy: The Data Says It Has Become a Global Shift. Healthsparkpr. 2026. Zdroj: <https://www.healthsparkpr.com/facial-ems-is-no-longer-a-fringe-therapy-the-data-says-it-has-become-a-global-shift/>
7. Hong, C. E., Park, H. R., & Kim, K. Y. (2017). Purchase Behavior, Usage Behavior and Satisfaction of Home Use Skin Care Devices during Home Shopping: Focus on 20's Age Women. Journal of the Korean Society of Cosmetology, 23(5), 1059-1070. Zdroj: <https://www.e-jkc.org/m/journal/view.php?number=2950>

8. Yu, Jin-Ho & Hwang, Tae-Yeon & Jeong, Jin-Gyu. (2009). The Influence of Medium Frequency Currents According to Frequency on Sensory Threshold. *Journal of the Korean Academy of Clinical Electrophysiology*. 7. 23-27. 10.5627/KACE.2009.7.1.023. Zdroj: https://www.researchgate.net/publication/272882724_The_Influence_of_Medium_Frequency_Currents_According_to_Frequency_on_Sensory_Threshold
9. Shin D, Park M, Choi Y, et al. Effects of medium frequency EMS (electrical muscle stimulation) At-home beauty device on the facial skin improvement. *J Korean Soc Cosmetol*. 2022;28:1031-1041. Zdroj: https://www.researchgate.net/publication/365128194_Effects_of_Medium_Frequency_EMS_Electrical_Muscle_Stimulation_At-home_Beauty_Device_on_the_Facial_Skin_Improvement
10. Park, J. S., Lee, J. W., & Lee, J. Y. (2011). The Effects of IFC and Ultrasound on the ROM and Pain in Patients with Chronic Back Pain. *Journal of the Korean Academy of Clinical Electrophysiology*, 9(1), 23-28. Zdroj: https://www.researchgate.net/publication/314585292_The_Effects_of_IFC_and_Ultrasound_on_the_ROM_and_Pain_in_Patients_with_Chronic_Back_Pain
11. Kim, C. Y. (2021). Effects of Interfering Current Stimulation on Vastus Medialis Oblique and Vastus Lateralis Activity and Ratio during Squat Exercise. *Journal of The Korean Society of Integrative Medicine*, 9(4), 283-289. Zdroj: <https://koreascience.kr/article/JAKO202134234725531.view>
12. Buschmann, I. Rubio, M.B., Barile, A.C., Grinberg, D. Modulated EMS in Your Fifties: A Post-Perimenopause Protocol. PureliftLab. 2026. Zdroj: <https://pureliftlab.com/blogs/skincare/modulated-ems-fifties-post-perimenopause>
13. Heidland A, Fazeli G, Klassen A, Sebekova K, Hennemann H, Bahner U, Di Iorio B. Neuromuscular electrostimulation techniques: historical aspects and current possibilities in treatment of pain and muscle wasting. *Clin Nephrol*. 2013 Jan;79 Suppl 1:S12-23. Zdroj: <https://pubmed.ncbi.nlm.nih.gov/23249528/>
14. Xu X, Zhang H, Yan Y, Wang J, Guo L. Effects of electrical stimulation on skin surface. *Acta Mech Sin*. 2021;37(12):1843-1871. doi: 10.1007/s10409-020-01026-2. Epub 2021 Feb 6. PMID: 33584001; PMCID: PMC7866966. Zdroj: <https://pmc.ncbi.nlm.nih.gov/articles/PMC7866966/>
15. Kwak, S.J et al. Effect of At-home Beauty Device Equipped with Microcurrent and Electrical Muscle Stimulation Functions on the Elasticity and Lifting of Facial Skin. *Journal of the Korean Society of Cosmetology* 2023;29(6):1505-1514. DOI: <https://doi.org/10.52660/JKSC.2023.29.6.1505> Zdroj: https://www.researchgate.net/publication/377267143_Effect_of_At-home_Beauty_Device_Equipped_with_Microcurrent_and_Electrical_Muscle_Stimulation_Functions_on_the_Elasticity_and_Lifting_of_Facial_Skin
16. Shin, Da-Som & Park, Min-Jung & Choi, Ye-Ji & Hur, Yu-Kyung. (2022). Effects of Medium Frequency EMS (Electrical Muscle Stimulation) At-home Beauty Device on the Facial Skin Improvement. *Journal of the Korean Society of Cosmetology*. 28. 1031-1041. 10.52660/JKSC.2022.28.5.1031. Zdroj: https://www.researchgate.net/publication/365128194_Effects_of_Medium_Frequency_EMS_Electrical_Muscle_Stimulation_At-home_Beauty_Device_on_the_Facial_Skin_Improvement
17. Brian R. The physiology of neuromuscular electrical stimulation. *Pediatr. Phys. Therapy*. 1997;9:96–102. Zdroj: <https://www.semanticscholar.org/paper/The-Physiology-of-Neuromuscular-Electrical-Reed/ad205f02aa364ba292892897b1728a58ac8fdb7c>
18. Kavanagh, S., Newell, J., Hennessy, M. and Sadick, N. (2012), Use of a neuromuscular electrical stimulation device for facial muscle toning: a randomized, controlled trial. *J Cosmet Dermatol*, 11: 261-266. <https://doi.org/10.1111/jocd.12007> Zdroj: <https://onlinelibrary.wiley.com/doi/abs/10.1111/jocd.12007>
19. Kemmler W, Schliffka R, Mayhew JL, von Stengel S. Effects of whole-body electromyostimulation on resting metabolic rate, body composition, and maximum strength in postmenopausal women: the Training and ElectroStimulation Trial. *J Strength Cond Res*. 2010 Jul;24(7):1880-7. doi: 10.1519/JSC.0b013e3181ddaeec. PMID: 20555279. Zdroj: <https://pubmed.ncbi.nlm.nih.gov/20555279/>
20. Yeom, S. Lee, H., & Jeon, K. (2022). *Korean Journal of Sport Science*, 33(1), 19-32, <https://doi.org/10.24985/kjss.2022.33.1.19>. Zdroj: <https://kjss.sports.re.kr/v.33/1/19/145>
21. EMS Revolution® Team. EMS TRAINING News. EMS TRAINING course. INDICATIONS AND CONTRAINDICATIONS. 2000. Zdroj: <https://www.emsrevolution.com/en/indications-and->

contraindications/

Rybná 1066/15, Staré Město, 110 00 Praha 1

www.bloominds.cz